

Product Datasheet

Mouse IL23 ELISA Kit (orb1202587)

Description	Mouse Interleukin 23 (IL23) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9EQ14
Sensitivity	5.1pg/mL
Entrez	83430
Expiration Date	Please enquire.