

Product Datasheet

Bovine Postacrosomal Sheath WW Domain Binding Protein (PAWP) ELISA Kit (orb1202613)

Description	Bovine Postacrosomal Sheath WW Domain Binding Protein (PAWP) ELISA Kit
Reactivity	Bovine
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	A3KFF6
Sensitivity	0.059ng/mL
Entrez	523070
Expiration Date	Please enquire.