

Product Datasheet

Bovine Xanthine Dehydrogenase (XDH) ELISA Kit (orb1202617)

Description	Bovine Xanthine Dehydrogenase (XDH) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P80457
Sensitivity	6.3pg/mL
Entrez	280960
Expiration Date	Please enquire.