

## Product Datasheet

### Rat Semaphorin 3A (SEMA3A) ELISA Kit (orb1202643)

|                          |                                                                     |
|--------------------------|---------------------------------------------------------------------|
| <b>Description</b>       | Rat Semaphorin 3A (SEMA3A)<br>ELISA Kit                             |
| <b>Reactivity</b>        | Rat                                                                 |
| <b>Range</b>             | 0.156-10ng/mL                                                       |
| <b>Note</b>              | For research use only                                               |
| <b>Application notes</b> | sandwich                                                            |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates or other biological<br>fluids. |
| <b>Uniprot ID</b>        | <b>Q63548</b>                                                       |
| <b>Sensitivity</b>       | 0.057ng/mL                                                          |
| <b>Entrez</b>            | <b>29751</b>                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                     |