

Product Datasheet

Rat Potassium Inwardly Rectifying Channel Subfamily J, Member 5 (KCNJ5) ELISA Kit (orb1202656)

Description	Rat Potassium Inwardly Rectifying Channel Subfamily J, Member 5 (KCNJ5) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P48548
Sensitivity	0.055ng/mL
Entrez	29713
Expiration Date	Please enquire.