

## Product Datasheet

# Rat Gamma-Secretase Activating Protein (gSAP) ELISA Kit (orb1202659)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Rat Gamma-Secretase<br>Activating Protein (gSAP)<br>ELISA Kit |
| <b>Reactivity</b>        | Rat                                                           |
| <b>Range</b>             | 0.312-20ng/mL                                                 |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | sandwich                                                      |
| <b>Sample Types</b>      | tissue homogenates or other<br>biological fluids.             |
| <b>Uniprot ID</b>        | <b>D3ZPQ9</b>                                                 |
| <b>Sensitivity</b>       | 0.113ng/mL                                                    |
| <b>Entrez</b>            | <b>311984</b>                                                 |
| <b>Expiration Date</b>   | Please enquire.                                               |