

Product Datasheet

Rat Fibronectin Type III Domain Containing Protein 5 (FNDC5) ELISA Kit (orb1202661)

Description	Rat Fibronectin Type III Domain Containing Protein 5 (FNDC5) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8K3V5
Sensitivity	0.59ng/mL
Entrez	260327
Expiration Date	Please enquire.