

Product Datasheet

Human EGF Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit (orb1202665)

Description	Human EGF Like Domain Protein, Multiple 7 (EGFL7) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9UHF1
Sensitivity	14.8pg/mL
Entrez	51162
Expiration Date	Please enquire.