

Product Datasheet

Sheep Insulin Like Growth Factor 1 (IGF1) ELISA Kit (orb1202683)

Description	Sheep Insulin Like Growth Factor 1 (IGF1) ELISA Kit
Reactivity	Sheep
Range	0.234-15ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P10763
Sensitivity	0.115ng/mL
Entrez	443318
Expiration Date	Please enquire.