

Product Datasheet

Sheep Ghrelin (GHRL) ELISA Kit (orb1202684)

Description	Sheep Ghrelin (GHRL) ELISA Kit
Reactivity	Sheep
Range	0.156-10ng/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q30DT1
Sensitivity	0.045ng/mL
Entrez	442996
Expiration Date	Please enquire.