

Product Datasheet

Bovine Deoxyribonuclease I (DNASE1) ELISA Kit (orb1202725)

Description	Bovine Deoxyribonuclease I (DNASE1) ELISA Kit
Reactivity	Bovine
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P00639
Sensitivity	0.32ng/mL
Entrez	282217
Expiration Date	Please enquire.