

Product Datasheet

Porcine Bcl2 Associated X Protein (Bax) ELISA Kit (orb1202729)

Description	Porcine Bcl2 Associated X Protein (Bax) ELISA Kit
Reactivity	Porcine
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q2A135
Sensitivity	0.28ng/mL
Entrez	396633
Expiration Date	Please enquire.