

Product Datasheet

Mouse Wingless Type MMTV Integration Site Family, Member 10B (WNT10B) ELISA Kit (orb1202753)

Description	Mouse Wingless Type MMTV Integration Site Family, Member 10B (WNT10B) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P48614
Sensitivity	0.056ng/mL
Entrez	22410
Expiration Date	Please enquire.