

Product Datasheet

Mouse NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit (orb1202768)

Description	Mouse NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q64669
Sensitivity	0.127ng/mL
Entrez	18104
Expiration Date	Please enquire.