

Product Datasheet

Mouse Nei Endonuclease VIII Like Protein 3 (NEIL3) ELISA Kit (orb1202770)

Description	Mouse Nei Endonuclease VIII Like Protein 3 (NEIL3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8K203
Sensitivity	0.061ng/mL
Entrez	234258
Expiration Date	Please enquire.