

Product Datasheet

Mouse Leucine Rich Repeat Containing Protein 3C (LRRC3C) ELISA Kit (orb1202775)

Description	Mouse Leucine Rich Repeat Containing Protein 3C (LRRC3C) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Sensitivity	0.068ng/mL
Entrez	192852
Expiration Date	Please enquire.