

Product Datasheet

Mouse Fibronectin Type III Domain Containing Protein 5 (FNDC5) ELISA Kit (orb1202784)

Description	Mouse Fibronectin Type III Domain Containing Protein 5 (FNDC5) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q8K4Z2
Sensitivity	0.246ng/mL
Entrez	384061
Expiration Date	Please enquire.