

Product Datasheet

Mouse Egl Nine Homolog 1 (EGLN1) ELISA Kit (orb1202791)

Description	Mouse Egl Nine Homolog 1 (EGLN1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q91YE3
Sensitivity	0.108ng/mL
Entrez	112405
Expiration Date	Please enquire.