

Product Datasheet

Mouse Deoxyribonuclease I Like Protein 3 (DNASE1L3) ELISA Kit (orb1202793)

Description	Mouse Deoxyribonuclease I Like Protein 3 (DNASE1L3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O55070
Sensitivity	0.058ng/mL
Entrez	13421
Expiration Date	Please enquire.