

Product Datasheet

Equine Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) ELISA Kit (orb1202813)

Description	Equine Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) ELISA Kit
Reactivity	Equine
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O77717
Sensitivity	1.49ng/mL
Entrez	100034134
Expiration Date	Please enquire.