

Product Datasheet

Equine Matrix Metalloproteinase 13 (MMP13) ELISA Kit (orb1202821)

Description	Equine Matrix Metalloproteinase 13 (MMP13) ELISA Kit
Reactivity	Equine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O18927
Sensitivity	6.9pg/mL
Entrez	100009711
Expiration Date	Please enquire.