

Product Datasheet

Goat Adrenocorticotrophic Hormone (ACTH) ELISA Kit (orb1202834)

Description	Goat Adrenocorticotrophic Hormone (ACTH) ELISA Kit
Reactivity	Goat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	5.15pg/mL
Expiration Date	Please enquire.