

Product Datasheet

Human Zinc Finger Homeobox Protein 4 (ZFHX4) ELISA Kit (orb1202837)

Description	Human Zinc Finger Homeobox Protein 4 (ZFHX4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q86UP3
Sensitivity	0.115ng/mL
Entrez	79776
Expiration Date	Please enquire.