

Product Datasheet

Human Zinc Finger, AN1-Type Domain Protein 6 (ZFAND6) ELISA Kit (orb1202838)

Description	Human Zinc Finger, AN1-Type Domain Protein 6 (ZFAND6) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q6FIF0
Sensitivity	0.056ng/mL
Entrez	54469
Expiration Date	Please enquire.