

## Product Datasheet

### Human Zinc Finger, BED Type Containing Protein 3 (ZBED3) ELISA Kit (orb1202839)

|                          |                                                                          |
|--------------------------|--------------------------------------------------------------------------|
| <b>Description</b>       | Human Zinc Finger, BED Type<br>Containing Protein 3 (ZBED3)<br>ELISA Kit |
| <b>Reactivity</b>        | Human                                                                    |
| <b>Range</b>             | 0.312-20ng/mL                                                            |
| <b>Note</b>              | For research use only                                                    |
| <b>Application notes</b> | sandwich                                                                 |
| <b>Sample Types</b>      | tissue homogenates, cell<br>lysates or other biological<br>fluids.       |
| <b>Uniprot ID</b>        | <b>Q96IU2</b>                                                            |
| <b>Sensitivity</b>       | 0.128ng/mL                                                               |
| <b>Entrez</b>            | <b>84327</b>                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                          |