

Product Datasheet

Human WW Domain Containing E3 Ubiquitin Protein Ligase 2 (WWP2) ELISA Kit (orb1202840)

Description	Human WW Domain Containing E3 Ubiquitin Protein Ligase 2 (WWP2) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	000308
Sensitivity	0.27ng/mL
Entrez	11060
Expiration Date	Please enquire.