

Product Datasheet

Human Wingless Type MMTV Integration Site Family, Member 3 (WNT3) ELISA Kit (orb1202847)

Description	Human Wingless Type MMTV Integration Site Family, Member 3 (WNT3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P56703
Sensitivity	0.058ng/mL
Entrez	7473
Expiration Date	Please enquire.