

Product Datasheet

Human Wingless Type MMTV Integration Site Family, Member 3A (WNT3A) ELISA Kit (orb1202848)

Description	Human Wingless Type MMTV Integration Site Family, Member 3A (WNT3A) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P56704
Sensitivity	0.062ng/mL
Entrez	89780
Expiration Date	Please enquire.