

Product Datasheet

Human Wingless Type MMTV Integration Site Family, Member 2 (WNT2) ELISA Kit (orb1202849)

Description	Human Wingless Type MMTV Integration Site Family, Member 2 (WNT2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P09544
Sensitivity	0.062ng/mL
Entrez	7472
Expiration Date	Please enquire.