

Product Datasheet

Human Wingless Type MMTV Integration Site Family, Member 2B (WNT2B) ELISA Kit (orb1202850)

Description	Human Wingless Type MMTV Integration Site Family, Member 2B (WNT2B) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q93097
Sensitivity	0.059ng/mL
Entrez	7482
Expiration Date	Please enquire.