

## Product Datasheet

### Human Testican 2 (TIC2) ELISA Kit (orb1202871)

|                          |                                                                  |
|--------------------------|------------------------------------------------------------------|
| <b>Description</b>       | Human Testican 2 (TIC2) ELISA Kit                                |
| <b>Reactivity</b>        | Human                                                            |
| <b>Range</b>             | 0.312-20ng/mL                                                    |
| <b>Note</b>              | For research use only                                            |
| <b>Application notes</b> | sandwich                                                         |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates<br>or other biological fluids. |
| <b>Uniprot ID</b>        | <b>Q92563</b>                                                    |
| <b>Sensitivity</b>       | 0.132ng/mL                                                       |
| <b>Entrez</b>            | <b>9806</b>                                                      |
| <b>Expiration Date</b>   | Please enquire.                                                  |