

Product Datasheet

Human Reactive Oxygen Species Modulator 1 (ROMO1) ELISA Kit (orb1202903)

Description	Human Reactive Oxygen Species Modulator 1 (ROMO1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P60602
Sensitivity	0.053ng/mL
Entrez	140823
Expiration Date	Please enquire.