

Product Datasheet

Human Ras Related C3 Botulinum Toxin Substrate 1 (Rac1) ELISA Kit (orb1202914)

Description	Human Ras Related C3 Botulinum Toxin Substrate 1 (Rac1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P63000
Sensitivity	0.060ng/mL
Entrez	5879
Expiration Date	Please enquire.