

Product Datasheet

Human RAB5C, Member RAS Oncogene Family (RAB5C) ELISA Kit (orb1202915)

Description	Human RAB5C, Member RAS Oncogene Family (RAB5C) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P51148
Sensitivity	0.055ng/mL
Entrez	5878
Expiration Date	Please enquire.