

Product Datasheet

Human Proline/Serine Rich Coiled Coil protein 1 (PSRC1) ELISA Kit (orb1202919)

Description	Human Proline/Serine Rich Coiled Coil protein 1 (PSRC1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q6PGN9
Sensitivity	0.048ng/mL
Entrez	84722
Expiration Date	Please enquire.