

Product Datasheet

Human Phospholipase A2, Group IVD (PLA2G4D) ELISA Kit (orb1202927)

Description	Human Phospholipase A2, Group IVD (PLA2G4D) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q86XP0
Sensitivity	0.064ng/mL
Entrez	280748
Expiration Date	Please enquire.