

Product Datasheet

Human Peptidase Inhibitor 16 (PI16) ELISA Kit (orb1202933)

Description	Human Peptidase Inhibitor 16 (PI16) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q6UXB8
Sensitivity	0.34ng/mL
Entrez	221476
Expiration Date	Please enquire.