

Product Datasheet

Human Parkinson Disease Protein 7 (PARK7) ELISA Kit (orb1202942)

Description	Human Parkinson Disease Protein 7 (PARK7) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q99497
Sensitivity	0.64ng/mL
Entrez	11315
Expiration Date	Please enquire.