

Product Datasheet

Human Nuclear Factor, Erythroid Derived 2 Like Protein 2 (NFE2L2) ELISA Kit (orb1202954)

Description	Human Nuclear Factor, Erythroid Derived 2 Like Protein 2 (NFE2L2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q16236
Sensitivity	0.051ng/mL
Entrez	4780
Expiration Date	Please enquire.