

Product Datasheet

Human Leucine Rich Repeat Containing Protein 3C (LRRC3C) ELISA Kit (orb1202980)

Description	Human Leucine Rich Repeat Containing Protein 3C (LRRC3C) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	A6NJW4
Sensitivity	0.066ng/mL
Entrez	100505591
Expiration Date	Please enquire.