

Product Datasheet

Human Longevity Assurance Homolog 2 (LASS2) ELISA Kit (orb1202981)

Description	Human Longevity Assurance Homolog 2 (LASS2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96G23
Sensitivity	0.061ng/mL
Entrez	29956
Expiration Date	Please enquire.