

Product Datasheet

Human Potassium Inwardly Rectifying Channel Subfamily J, Member 10 (KCNJ10) ELISA Kit (orb1202990)

Description	Human Potassium Inwardly Rectifying Channel Subfamily J, Member 10 (KCNJ10) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P78508
Sensitivity	0.64ng/mL
Entrez	3766
Expiration Date	Please enquire.