

Product Datasheet

Human G Protein Coupled Receptor, Family C, Group 5, Member A (GPRC5A) ELISA Kit (orb1203007)

Description	Human G Protein Coupled Receptor, Family C, Group 5, Member A (GPRC5A) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8NFJ5
Sensitivity	0.116ng/mL
Entrez	9052
Expiration Date	Please enquire.