

Product Datasheet

Human Growth Arrest And DNA Damage Inducible Protein Alpha (GADD45a) ELISA Kit (orb1203014)

Description	Human Growth Arrest And DNA Damage Inducible Protein Alpha (GADD45a) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P24522
Sensitivity	0.064ng/mL
Entrez	1647
Expiration Date	Please enquire.