

Product Datasheet

Human ER Lipid Raft Associated Protein 2 (ERLIN2) ELISA Kit (orb1203024)

Description	Human ER Lipid Raft Associated Protein 2 (ERLIN2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	O94905
Sensitivity	0.115ng/mL
Entrez	11160
Expiration Date	Please enquire.