

Product Datasheet

Human Defensin Beta 131 (DEFb131) ELISA Kit (orb1203041)

Description	Human Defensin Beta 131 (DEFb131) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P59861
Sensitivity	5.7pg/mL
Entrez	644414
Expiration Date	Please enquire.