

Product Datasheet

Human Cytochrome P450 2C19 (CYP2C19) ELISA Kit (orb1203048)

Description	Human Cytochrome P450 2C19 (CYP2C19) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P33261
Sensitivity	0.056ng/mL
Entrez	1557
Expiration Date	Please enquire.