

Product Datasheet

Human CUB And Zona Pellucida Like Domains Protein 1 (CUZD1) ELISA Kit (orb1203049)

Description	Human CUB And Zona Pellucida Like Domains Protein 1 (CUZD1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q86UP6
Sensitivity	0.114ng/mL
Entrez	50624
Expiration Date	Please enquire.