

## Product Datasheet

### Human Collagen Triple Helix Repeat Containing Protein 1 (CTHRC1) ELISA Kit (orb1203050)

|                          |                                                                            |
|--------------------------|----------------------------------------------------------------------------|
| <b>Description</b>       | Human Collagen Triple Helix Repeat Containing Protein 1 (CTHRC1) ELISA Kit |
| <b>Reactivity</b>        | Human                                                                      |
| <b>Range</b>             | 1.56-100ng/mL                                                              |
| <b>Note</b>              | For research use only                                                      |
| <b>Application notes</b> | sandwich                                                                   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates or other biological fluids.              |
| <b>Uniprot ID</b>        | <b>Q96CG8</b>                                                              |
| <b>Sensitivity</b>       | 0.58ng/mL                                                                  |
| <b>Entrez</b>            | <b>115908</b>                                                              |
| <b>Expiration Date</b>   | Please enquire.                                                            |