

Product Datasheet

Human C-Type Lectin Domain Family 11, Member A (CLEC11A) ELISA Kit (orb1203060)

Description	Human C-Type Lectin Domain Family 11, Member A (CLEC11A) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9Y240
Sensitivity	0.51ng/mL
Entrez	6320
Expiration Date	Please enquire.